



... GREEN AIR™

UNIQUE BREATHABLE GLOVES
FOR MAXIMUM COMFORT



Compostable glove: The one and only



..... High degree of breathability & personal safety

- Water vapour passes 7 more through than classic gloves for a sweat-free experience. GreenGlove Air comfort: 221g/m²/day vs 31g/m²/day for PE gloves.
- Safe material: biocompatible and non-cytotoxic for less irritations to the skin.

..... Unique characteristics suitable for gentle and natural calving, foaling and farrowing

- Soft and strong: Highly stretchable and resistant to puncture.
- Tactile: For better handling, feel and touch of the professional.
- Natural: An earth-based product to greet your animal into this world.



Compostable and biodegradable

Sustainable development, an IMV Group commitment



The IMV Technologies Group is committed to sustainable development. This glove is the fruit of a constant focus and attention on the environment whereby we search for raw materials or production methods that have a minimal or no impact on the environment. Beyond products, we also motivate our employees in such a way that the environmental aspects become a natural part of our workday.

- An environmental alternative to polyethylene-based plastics, GreenGlove™ is the first completely biodegradable and compostable long glove ever produced.
- Reduced greenhouse gas emissions, consumption of energy and non-renewable resources during its manufacturing, thanks to the use of raw materials of agricultural origin.

Composted in 6 months.

Ref. 026066 - Air Comfort GreenGlove

Glove with tear-off strip

One-size (L)

Length 920 mm (36.5") Thickness 21 microns (0.85 mil")

Box of 100 gloves - Parcel of 10 x 100 gloves

